


# Lx Video Capture

The **Lx Video Capture** program assists work in the clinic by enabling the capture, display, storage, retrieval, editing and printing of video images on a PC, desktop or laptop.

High quality digital pictures, storage on hard disk or other PC media, instantaneous viewing, copy/paste into reports and printing on standard printers are all immediately available.



The **Lx Video Capture** system can be used with any continuous or stroboscopic light source, with rigid or flexible endoscopes and with a range of cameras.

## Applications

- Video stroboscopy
- Nasoendoscopy
- Otoscopy

**Lx**  
Laryngograph

ISO9001:2000  
ISO13485:2003  
FM73968

## Features

Brings the power of the PC to all situations where direct observation or VCRs are now used.

Works with any combination of light source, camera and endoscope.

Captures images onto a desktop PC or a laptop for portability.

High quality video images are displayed on the PC monitor and captured to hard disk.

Slow motion animation of vocal fold cycles when used with a strobe light source.

Simultaneous capture and display of EGG, acoustic or other signals.

Footswitch control to start/stop recording, mark sections of interest, grab images and start/pause/stop playback.

Continuous or field-by-field viewing available immediately after recording ends.

Easy editing functions for production of concise image sequences.

Side-by-side replay of multiple examinations for evidence based outcomes

Easy archiving onto DVD.

**Laryngograph Ltd.** • 1 Foundry Mews, London, NW1 2PR  
Tel: +44 (0)20 7387 7793 • Fax: +44 (0)20 7383 2039  
E-mail: lx@laryngograph.com • www.laryngograph.com